

## Visualizing Biotinylated Dextran Amine (BDA) Protocol Chromogenic or Fluorescent

### Bleach Solution

(50% methanol & 1% hydrogen peroxide in PBS) 18 mL

- 30% Hydrogen peroxide ( $H_2O_2$ ) 700  $\mu$ L
- Methanol ( $CH_3OH$ ) 9 mL
- 1x PBS 9 mL

### Streptavidin Solution (per tray, ~20 mL)

0.4% PBS-Tx 20 mL

Streptavidin (fluorescent)

or

Streptavidin – HRP (chromogenic) \_\_\_\_  $\mu$ L

Streptavidin type	Storage [ ] (%)	Dilution factor	Volume of reagent ( $\mu$ L)	Volume of 0.4% PBS-Tx (mL)
_____	_____	1 : _____	_____	20

e.g. in order to get a dilution of 1:250 for streptavidin with 50% concentration:

- Multiply the dilution factor (1/250) with the volume of PBS-Tx (20)
- Convert the resulting volume from mL to  $\mu$ L by multiplying it by 1000
- Correct for the storage concentration by multiplying the volume by 2 (because 100/50)

### Protocol

*For each step: Make ~20 mL of solution per tray. Fill the wells up to around half-full so that the sections are fully submerged without being able to float out over the top. Make sure all sections are free-floating and not stuck to the sides. Put the tray on the rotator. Dry sections off between each wash by dabbing the tray on absorbent bench liner.*

1. Wash sections in 1x PBS for 5x5 minutes.
2. Bleach sections for 20 minutes.
3. Wash sections in 1x PBS for 5x5 minutes.
4. Incubate sections in streptavidin solution for 1 hour at room temp.
  - a. For fluorescent injections: from this point on, cover sections with tinfoil to reduce light exposure to fluorescent-stained tissue.
5. Wash sections in 1x PBS for 5x5 minutes.
6. Store the sections at 4°C in either 1x PBS (for short-term storage, <24 hrs) or 1x PBS Azide (for long-term storage, >24 hrs).

## DAB Staining Protocol

### Chromogenic

#### Diaminobenzidine Peroxidase Reaction Medium (for **BROWN** reaction product)

	50 mL	100 mL
– 30% Hydrogen peroxide ( $H_2O_2$ )	7.5 $\mu$ l	15 $\mu$ l
– DAB (Sigma)	12.5 mg	25 mg
– 1x PBS	50 mL	100 mL

#### Diaminobenzidine peroxidase reaction medium (for **BLACK** reaction product)

	50 mL	100 mL
– 0.5% Cobalt(II) chloride ( $CoCl_2$ )	1.5 mL	3 mL
– 30% Hydrogen peroxide ( $H_2O_2$ )	7.5 $\mu$ l	15 $\mu$ l
– DAB (Sigma)	12.5 mg	25 mg
– 1x PBS	48.5 mL	97 mL

#### Protocol

1. Prepare a waterbath containing bleach to neutralize spills and to deactivate DAB after experiment is finished.
2. Prepare diaminobenzidine peroxidase reaction medium.
  - a. BDA is detected with a black DAB reaction while CTB is detected with a brown DAB reaction.
3. Incubate the sections in the reaction medium until the reaction is considered complete – this can be from 10 seconds to 2 minutes.
4. Wash sections in 1x PBS for 3x10 minutes.