

## OCT-Embedding and Flash-Freezing Protocol

### Solutions needed

<i>Dry ice</i>	<i>2-3 lbs</i>
<i>Optimal Cutting Temperature (OCT) compound</i>	<i>as needed</i>

### Materials needed

Falcon tube	Red plastic tub & lid
Silicon ice cube tray/mold	Large polystyrene cooler
Metal forceps (sharp or blunt)	Paper towel

### Protocol

1. Bring the red plastic tub and lid to chem stores to pick up dry ice.
  - a. Chem stores is open from 9-11:30am and 12:30-4pm Monday through Friday.
2. Empty dry ice from red tub into large polystyrene cooler then shake the cooler to generate a flat surface for the silicone ice cube tray/mold to sit atop and remain level and undistorted.
3. Pour a small amount of OCT to one silicone ice cube mold, only use enough to fully coat the bottom of the mold.
4. Transfer the silicone tray/mold into the polystyrene cooler.
  - a. Make sure the tray and the liquid within it are sitting level atop the dry ice.
  - b. Replace the lid and wait ~5-10 minutes for OCT to freeze (OCT will turn white once it is frozen).
5. Weigh the brain just prior to starting step 6.

**\*\* The remaining steps must be performed quickly to prevent the OCT from thawing. \*\***

6. Remove the silicone tray/mold from the polystyrene cooler.
  - a. Add OCT to mold containing frozen base layer, fill to approximately 2/3 full and avoid producing air bubbles.
  - b. Use a metal forceps to remove and/or move any bubbles away to the side.
7. Place the brain into the OCT.
  - a. Use metal forceps to move the brain to the center of the mold (diagonal with the caudal side up) See photo 1 for how to orient the brain. Remove and/or move any bubbles away to the side.
  - b. Pour OCT on top of brain, ensuring that it is fully submerged and covered in OCT.
  - c. Take note of the orientation of the brain and which corner the dorsal side is on.
8. Put the ice cube tray back into the polystyrene cooler for ~5-10 minutes, until block is completely frozen.
9. Remove the tray. Using a sharpie, mark the corner where the dorsal side of the brain is facing.
10. Remove the OCT-embedded brain from the silicone ice cube tray/mold by everting the mold into an empty falcon tube.
11. Place the falcon tube containing the OCT-embedded brain into the -80°C freezer.

**Photo 1.**

Top view – caudal side shown

